



FEMUR

KEEL AND SINGLE FEMORAL POST

- Provides restoration of femoral condylar surface anatomy.
- Improves implant strength, stability and resistance to rotation through a mechanical lock.
- Peripheral rim to minimize cement extrusion



TIBIA

ALL-POLY INLAY TIBIA

- Provides restoration of tibial condylar surface anatomy.
- Component captured with optimized peripheral cement technique.

THE RENEW™ KNEE PROVIDES THE FOLLOWING:

- Resurfacing unicondylar knee platform designed for outpatient setting with same day discharge.
- Minimal bone resurfacing without opposite tibial-femoral or patellofemoral involvement.
- Avoids quad interruption, patellar eversion/subluxation and ligament sacrifice.
- Maximizes OR and surgical efficiencies.
- Reproducible surgical technique.
- Faster recovery and return to activities of daily living.
- Yields expected and proven results.

When appropriately diagnosed, resurfacing unicondylar knee arthroplasty (rUKA) avoids the finality of total knee arthroplasty.

TIBIAL COMPONENTS

SIZE	PN	DESCRIPTION	AP	HEIGHT
SMALL	311326	Partial All Poly Tibial Component, Small, 32mm, 6.5mm	32	6.5
	311327	Partial All Poly Tibial Component, Small, 32mm, 7.5mm		7.5
	311328	Partial All Poly Tibial Component, Small, 32mm, 8.5mm		8.5
SMALL-MEDIUM	311346	Partial All Poly Tibial Component, Sm-Med, 34mm, 6.5mm	34	6.5
	311347	Partial All Poly Tibial Component, Sm-Med, 34mm, 7.5mm		7.5
	311348	Partial All Poly Tibial Component, Sm-Med, 34mm, 8.5mm		8.5
MEDIUM	311366	Partial All Poly Tibial Component, Medium, 36mm, 6.5mm	36	6.5
	311367	Partial All Poly Tibial Component, Medium, 36mm, 7.5mm		7.5
	311368	Partial All Poly Tibial Component, Medium, 36mm, 8.5mm		8.5
MEDIUM-LARGE	311386	Partial All Poly Tibial Component, Med-Lrg, 38mm, 6.5mm	38	6.5
	311387	Partial All Poly Tibial Component, Med-Lrg, 38mm, 7.5mm		7.5
	311388	Partial All Poly Tibial Component, Med-Lrg, 38mm, 8.5mm		8.5
LARGE	311406	Partial All Poly Tibial Component, Large, 40mm, 6.5mm	40	6.5
	311407	Partial All Poly Tibial Component, Large, 40mm, 7.5mm		7.5
	311408	Partial All Poly Tibial Component, Large, 40mm, 8.5mm		8.5
EXTRA-LARGE	311426	Partial All Poly Tibial Component, Xlarge, 42mm, 6.5mm	42	6.5
	311627	Partial All Poly Tibial Component, Xlarge, 42mm, 7.5mm		7.5
	311428	Partial All Poly Tibial Component, Xlarge, 42mm, 8.5mm		8.5

FEMORAL COMPONENTS

SIZE	PN	DESCRIPTION	AP	ML	POST LENGTH
EXTRA SMALL	313421	Partial Femoral Component, Left, XSML, SZ 42x15	42	15	11
SMALL	313451	Partial Femoral Component, Left, SML, SZ 45x16	45	16	11
SMALL MEDIUM	313471	Partial Femoral Component, Left, SMED, SZ 47x17	47	17	11
MEDIUM	313501	Partial Femoral Component, Left, MED, SZ 50x18	50	18	11
LARGE	313541	Partial Femoral Component, Left, LRG, SZ 54x19	54	19	11
EXTRA SMALL	314421	Partial Femoral Component, Right, XSML, SZ 42x15	42	15	11
SMALL	314451	Partial Femoral Component, Right, SML, SZ 45x16	45	16	11
SMALL MEDIUM	314471	Partial Femoral Component, Right, SMED, SZ 47x17	47	17	11
MEDIUM	314501	Partial Femoral Component, Right, MED, SZ 50x18	50	18	11
LARGE	314541	Partial Femoral Component, Right, LRG, SZ 54x19	54	19	11

Freedom Total Knee® System (Renew Resurfacing Unicondylar)



Maxx Orthopedics, Inc.
2460 General Armistead Ave, Ste 100
Norristown, PA 19403 USA



AJW Consulting GmbH
Breite Straße 3
40213 Düsseldorf, Germany
Phone: +49(0) 211 54059 6030
Email: info@ajwtech.com



Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the Freedom Total Knee® System.

CAUTION: THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

WARNING: THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY

Freedom Knee® is manufactured by Maxx Orthopedics, Inc. FREEDOM, FREEDOM KNEE and RENEW KNEE are Registered Trademarks of Maxx Orthopedics, Inc.

© 2022 Maxx Orthopedics. All rights reserved. Updated February 2022

Maxx Orthopedics Brochure No. MXO-MP00048-R01